D

LAVANOX-Serag® Wundspüllösung 250ml

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>99% Aqua purificata, Natriumhypochlorit/hypochlorige Säure

GB

LAVANOX-Serag® Wound Irrigation Solution 250ml

Instructions for use

Intended purpose

For cleansing, moisturing and infection prophylaxis of acute and chronic wounds.

Product features

Wound cleansing and moistening with LAVANOX-Serag® wound irrigation solution helps to remove biofilm and fibrinoid debris. This promotes wound healing and reduces the risk of infection. The LAVANOX-Serag® wound irrigation solution also supports a low-pain removal of the wound dressing. Wound odours are eliminated by using LAVANOX-Serag® wound irrigation solution. The application of LAVANOX-Serag® wound irrigation solution is well tolerated and virtually painless. LAVANOX-Serag® wound irrigation solution is suitable for wounds on skin and mucous membranes. There are no restrictions on the frequency and duration of use. Residues of LAVANOX-Serag® wound irrigation solution may remain in the wound. As a multi-dose container, LAVANOX-Serag® wound irrigation solution is preserved with NaOCl/HOCl (<0.08%). This means that the opened single product can be used for three months. A slight chlorine odour is a sign of the activity of the preservative system.

Instructions for use

A wound should first be cleaned of

dirt, blood and coatings. For breastfeeding and in children under two years of age, consult your doctor or pharmacist beforehand.

Intended users
 Health care professionals and laypersons. LAVANOX-Serag® wound irrigation solution can be used by medical laypersons if they have been instructed by medical professionals. Children up to the age of 12 should be supervised by an adult.

Medical indications
 Acute and chronic wounds.


Limitations of use and contra-indications
 LAVANOX-Serag® wound irrigation solution should not be used in the case of known allergies to any of the ingredients. If you notice any signs of an allergic reaction to the ingredients, discontinue use of the product. Do not swallow.


Warnings, precautions, measures to be taken
 Avoid contact with eyes. In case of eye contact, rinse with plenty of water. Despite the generally very good tolerability, there is a residual risk of a discomforting feeling during application due to the individual wound situation. Allergic reactions cannot be ruled out. Changing the dressing too often (overdosing) can disrupt the wound healing process. The treatment should therefore be monitored by a healthcare professional. If there is redness, swelling or pain in the wound area or if the condition of the wound does not improve within a few days, a doctor should be consulted. Only use the product up to the stated expiry date and within 3 months after opening. Please report any potentially serious incident related to the product immediately to Pronotmed (vigilanz@pronotmed.de) and the competent authority in your country.

Patient target group(s)
 LAVANOX-Serag® wound irrigation solution can be used without restriction with regard to the age and gender of the patient. For use during pregnancy, while

Safe disposal
 No special measures are required for safe disposal. Please observe the local disposal guidelines.

Composition
 >99% Aqua purificata, sodium hypochlorite/hypochlorous acid

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GmbH & Co. KG
Pronotmed GmbH