D

LAVANOX-Serag® Wundgel 50ml

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Instructions for use

Intended purpose

For cleansing, long-term moistu-
rising and infection prophylaxis of acute and chronic wounds.

Product features

LAVANOX-Serag® wound gel is a ready-to-use hydrogel for painless application to acute and chronic wounds. It forms a protective film and keeps the wound moisturised. It thus forms a microenvironment that promotes wound healing and supports the natural healing and cleansing process in the wound. The risk of infection is reduced. Wound odours are eliminated by the use of LAVANOX-Serag® wound gel. The application of LAVANOX-Serag® wound gel is well tolerated and virtually painless. LAVANOX-Serag® wound gel is suitable for wounds on skin and mucous membranes. There are no restrictions on the frequency and duration of use. Residues of LAVANOX-Serag® wound gel may remain in the wound. As a multi-dose container, LAVANOX-Serag® wound gel is preserved with NaOCl/HOCl (<0.06%). This means that the opened single product can be used for three months. A slight chlorine odour is a sign of the activity of the preservative system.

Instructions for use

Shake well before use! Before the first application and when changing the dressing, we recommend using a wound spray or wound irrigation solution to remove exudate, fibrinoid deposits and other impurities. Products from the LAVANOX-Serag® product line are suitable for this purpose. Generously cover the wound surface with the LAVANOX-Serag® wound gel. Dress

the wound for protection. No undesirable side effects are known. The frequency of dressing changes depends on the nature of the wound and must be adapted to the respective wound situation. LAVANOX-Serag® wound gel can be used in combination with conventional dressings, plasters and other bandages. The simultaneous use of LAVANOX-Serag® wound gel with pre-soaked wound dressings, antiseptics and other wound care products can weaken their effect. Interactions are possible with chlorhexidine and polyurethane. Dressing materials should therefore be compatible with oxidising agents, especially hypochlorite. However, delayed application is possible. Please also refer to the instructions for use of the respective product.

Patient target group(s)

LAVANOX-Serag® wound gel can be used without restriction with regard to the age and gender of the patient. For use during pregnancy, while breastfeeding and in children under two years of age, consult your doctor or pharmacist beforehand.

Intended users

Health care professionals and laypersons. LAVANOX-Serag® wound gel can be used by medical laypersons if they have been instructed by medical professionals. Children up to the age of 12 should be supervised by an adult.

Medical indications

Acute and chronic wounds.

Limitations of use and contra-indications

LAVANOX-Serag® wound gel should not be used in the case of known allergies to any of the ingredients. If you notice any signs of an allergic reaction to the ingredients, discontinue use of the product. Do not swallow. Do not use gel that has been frozen. Do not use LAVANOX-Serag® wound gel on bleeding or heavily oozing wounds.

Warnings, precautions, measures to be taken


Avoid contact with eyes. In case of eye contact, rinse with plenty of water. Despite the generally very good tolerability, there is a residual risk of a discomforting feeling during application due to the individual wound situation. Allergic reactions cannot be ruled out. Changing the dressing too often (overdosing) can disrupt the wound healing process. The treatment should therefore be monitored by a healthcare professional. If there is redness, swelling or pain in the wound area or if the condition of the wound does not improve within a few days, a doctor should be consulted. Only use the product up to the stated expiry date and within 3 months after opening. Please report any potentially serious incident related to the product immediately to Promotomed (vigilanz@promotomed.de) and the competent authority in your country.




Safe disposal

No special measures are required for safe disposal. Please observe the local disposal guidelines.

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GmbH & Co. KG
Promotomed GmbH