## Artwork Info:

Product Description: **Iodosorb Dressing IFU**, Component: **300 x 600 mm**

Global Dimensions: **Panel 50 x 300mm**

Product Code: **66001290/66001291** Colours: **Black**  
**66001292/66001293**

Artwork Ref/BSC: **46275 / 0350209**

Prev Artwork Ref/BSC: **21199**

Date: **30 Jan 2025**

Revision: **02** Created by: **PK**

<b>(EN)</b>	<b>SmithNephew</b>
<b>IODOSORB® DRESSING</b>	<b>IODOSORB® DRESSING</b>
Cadexomer dressing with iodine	
<b>Description</b>	
IODOSORB® Dressings are sterile dark brown paste dressings with a gauze backing on both sides. The paste consists of cadexomer (modified starch microbeads), polyethylene glycol and iodine. The cadexomer beads are biodegradable. IODOSORB Dressings will absorb excess exudate and slough from the wound bed, and kill bacteria at the wound surface. In doing this it transforms into a soft moist gel.	
<b>Indications</b>	
IODOSORB Dressings are indicated for the topical treatment of chronic exuding wounds such as venous leg ulcers, diabetic foot ulcers and pressure injuries. The ulcers must be free of infection and IODOSORB Dressings can be used under compression therapy. IODOSORB Dressings may be used on infected wounds. Where the product is used on infected wounds, it is not intended for use as a dressing and treated as per local clinical protocol.	
<b>Contraindications</b>	
IODOSORB Dressings should not be used on dry necrotic tissue, or on patients with a known sensitivity to iodine or any of the other components. Do not use IODOSORB Dressings on children, pregnant or lactating women or patients with thyroid disorders (e.g. Hashimoto's thyroiditis or non-toxic nodular goitre).	
<b>Warnings</b>	
There is a potential risk of interaction with lithium, resulting in an increased possibility of hypothyroidism. Do not use IODOSORB Dressings concomitantly with mercurial antiseptics or taurine. If the patient is taking a lithium therapy functional tests, they should inform the healthcare professional.	
<b>Precautions for use</b>	
The maximum dose of the product must be observed. Particular care must be taken when using in patients with renal impairment or patients who may develop renal impairment, especially when treating large wounds or during prolonged use. IODOSORB Dressings are single use only and should not be re-used as cross-contamination or infection may occur and the dressing may be ineffective. Do not use IODOSORB Dressings in the vicinity of the eyes, nose or mouth. The product is for external use only. IODOSORB Dressings may cause a slight transient pain in the first hour after treatment. This is a sign that the product is beginning to clean the wound. Occasionally, IODOSORB Dressings may cause the skin around the wound edges to swell or redden. This will usually pass. If these symptoms persist or if the patient experiences any other symptoms, they should contact the healthcare professional. Contact with the skin around the wound edges/infect skin should be minimised. It is possible for an adherent crust to form if the IODOSORB Dressings are allowed to dry out. A single application of IODOSORB Dressings should not exceed 50g (equivalent to 5 x 10g sachets) and not more than 150g (equivalent to 15 x 10g sachets) applied in one treatment session. The duration of treatment should not exceed 3 weeks. The product is for healthcare professional use only.	
<b>Instructions for use</b>	
<b>Applying IODOSORB Dressing</b>	
1. Wash your hands thoroughly before and after use. Using gloved hands, clean the wound and the surrounding area with either a gentle stream of water or saline. DO NOT DRY the wound surface. 2. Remove the carrier gauge on both sides of the past. 3. Cut or mould to the shape of the wound. Cover the wound completely with an appropriate secondary dressing. Apply compression bandaging if appropriate. Any remaining IODOSORB Dressing should be discarded due to loss of sterility.	
<b>Changing IODOSORB Dressings</b>	
1. IODOSORB Dressings should be changed when it has become saturated with exudate and all the iodine has been released. This is indicated by loss of colour, usually two to three times a week. If the wound is discharging heavily, daily changes may be needed. If IODOSORB Dressing does not become saturated by exudate, the maximum time between changes is 3 days. 2. If necessary, soak the dressing for a few minutes, then remove. 3. Gently remove the IODOSORB Dressing using a stream of water or saline. 4. Gently blot any excess fluid, leaving the wound surface slightly moist, before re-applying IODOSORB Dressing. 5. Dispose of according to local clinical protocols, or throw away hygienically after use. If during the use of this device you believe that a serious incident has occurred, please report it to the manufacturer and your competent authority. Complaints and enquiries should be sent to: A copy of the summary of safety and clinical performance (SSCP) is available via https://ec.europa.eu/tools/eudamed	
<b>Product Availability</b>	
66001290 5x5g	
66001291 3x10g	
66001292 5x 10 g	
66001293 2x17 g	
<small>*Trade Mark of Smith-Nephew. Date of publication 01/2025. Revision 01. This product may be covered by one or more US patents. www.smith-nephew.com/patents</small>	

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