ACTICOAT[◊] FLEX 3

Silver-coated antimicrobial barrier dressing

Description

ACTICOAT[◊] Flex 3 Antimicrobial Barrier Dressings consist of a flexible, low adherent polyester layer coated with nanocrystalline silver. ACTICOAT Flex 3 Dressings are highly conformable dressings that follow the body contours to maintain contact with the wound surface. The dressings are low adherent, which helps to minimise wound trauma at dressing changes.

Mode of Action

ACTICOAT Flex 3 Dressings provide an effective antimicrobial barrier to microbial contamination protecting the wound from invasive pathogenic micro-organisms, thus assisting fast healing. Nanocrystalline silver is also effective against micro-organisms present in the wound, helping to reduce the risk of infection. The antimicrobial barrier properties of ACTICOAT Flex 3 Dressings remain effective for the wear time of the dressing.

Indications

ACTICOAT Flex 3 Dressings are antimicrobial barrier dressings indicated for use over partial and full thickness wounds such as burns, surgical sites, pressure ulcers, venous ulcers and diabetic ulcers. ACTICOAT Flex 3 Dressings may be used on infected wounds. Where the product is used on infected wounds, the infection should be inspected and treated as per local clinical protocol. ACTICOAT Flex 3 Dressings may be used as a wound contact layer in combination with Negative Pressure Wound Therapy (NPWT) for a period of up to 3 days.

Contraindications

- Do not use on patients with a known sensitivity to silver.
- Do not use on patients during MRI (Magnetic Resonance Imaging) examination.
- Prior to administering radiation therapy, remove ACTICOAT Flex 3 Dressing. A new dressing can be applied following each treatment.

Precautions

- For external use only, for example do not apply the dressing to exposed organs.

contamination or infection may occur and the dressings may be ineffective. In addition, the antimicrobial properties of the product will be compromised. Opening the dressing pack compromises the sterile barrier. Therefore, any unused dressings should not be retained for application at a later date.

ACTICOAT Flex 3 Dressings are for healthcare professional use only.

Instructions for use

- Where required, wound cleansing should be performed according to local clinical protocol. DO NOT use saline, as this may compromise the antimicrobial activity.
- Remove the ACTICOAT Flex 3 Dressing from the pack using a clean technique.
- Cut the dressing to shape as necessary. When used with compression therapy, the dressing must be cut to the size of the wound.
- Moisten the dressing with sterile water according to local clinical protocol, fully submerging the dressing then squeezing gently (DO NOT use saline). Allow excess water to drain prior to application. Alternatively a thin layer of a hydrogel product e.g. INTRASITE[◊] Gel can be applied directly to the wound and/or the dressing.
- When in use with NPWT, the wound filler material should be moistened and changed according to manufacturer's instructions. When using ACTICOAT Flex 3 Dressings under NPWT, the use of hydrogels e.g. INTRASITE is not recommended.
- Without stretching the dressing, apply ACTICOAT Flex 3 Dressing either way up, to the wound surface ensuring there are no creases. The dressing should be applied with the direction of stretch running along the limb to allow movement.
- ACTICOAT Flex 3 Dressings may be used to pack wounds.
- Secure the ACTICOAT Flex 3 dressing in place using an appropriate secondary dressing or compression therapy that will maintain a moist wound environment. In the case of highly exuding wounds, an absorbent

- ACTICOAT Flex 3 Dressings are not compatible with oil based products such as petrolatum.
 - ACTICOAT Flex 3 Dressings may not be compatible with topical antimicrobials.
 - Avoid contact with electrodes or conductive gels during electronic measurements e.g. EEG and ECG.
 - There may be transient pain experienced on application of ACTICOAT Flex 3 Dressings. Should continuous pain be experienced after application, remove the dressing. To alleviate this pain, apply a new dressing and use a hydrogel product as described below. Should pain persist, remove the dressing and discontinue use.
 - ACTICOAT Flex 3 Dressings should only be used on premature infants (less than 37 weeks gestation) when the clinical benefit outweighs potential risks. No clinical data is available in this age group and only limited data is available for use in neonates.
 - Clinicians/Healthcare Professionals should be aware that there is very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
 - Avoid use of circumferential dressing application when using ACTICOAT FLEX 3 10cm x 120cm / 4in. x 48in dressings. Instead consider using multiple pieces to minimize risk of decreased distal circulation. Remove if changes in circulation are detected.
 - ACTICOAT Flex 3 Dressings are single use only and should not be re-used as cross secondary dressing is appropriate. ACTICOAT Flex 3 Dressing allows transport of the exudate directly into the secondary dressing.
 - Keep ACTICOAT Flex 3 Dressing moist, but not so wet that tissue maceration occurs.
 - Change the dressing depending on the amount of exudate present and the condition of the wound. Dressings may be left in place for up to 3 days. Dispose of according to local clinical protocols, or throw away hygienically after use.
 - When using ACTICOAT Flex 3 Dressings under compression, regularly check the entire dressing system and change only if strike through occurs. An additional non-occlusive absorbent dressing may be required if the level of exudate is high.
 - Should the dressing dry or adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.
- Note:** ACTICOAT Flex 3 Dressings may cause transient discolouration of the surrounding skin.
- Note:** The colour of the dressing may vary. This does not affect the performance of the dressing.
- If during the use of this device you believe that a serious incident has occurred, please report it to the manufacturer and your competent authority.
complaints@smith-nephew.com
 A copy of the summary of safety and clinical performance (SSCP) is available via <https://ec.europa.eu/tools/eudamed>

Product Availability

Code	Size	Availability
66800396	5cm x 5cm / 2in. x 2in.	Carton of 5
66800398	10cm x 10cm / 4in. x 4in.	Carton of 5
66800399	10cm x 10cm / 4in. x 4in.	Carton of 12
66800409	10cm x 20cm / 4in. x 8in.	Carton of 12
66800419	20cm x 40cm / 8in. x 16in.	Carton of 6
66800432	40cm x 40cm / 16in. x 16in.	Carton of 6
66800435	10cm x 120cm / 4in. x 48in.	Carton of 6
66801290	4cm x 15cm / 1 ⁹ / ₁₆ in. x 6in.	Carton of 5
66801291	4cm x 25cm / 1 ⁹ / ₁₆ in. x 10in.	Carton of 5
66801292	4cm x 35cm / 1 ⁹ / ₁₆ in. x 13 ³ / ₄ in.	Carton of 5
66801293	20cm x 20cm / 8in. x 8in.	Carton of 5

°Trade Mark of Smith+Nephew. ©Smith+Nephew 2024 Date of issue 04/2024.

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°Trade Mark of Smith+Nephew ©Smith+Nephew 2024
Smith & Nephew Medical Limited, 101 Hessle Road, Hull HU3 2BN England.
Smith & Nephew Operations B.V., Bloemlaan 2, 2132 NP Hoofddorp, Netherlands