D

LAVANOX-Serag® Wundspray 75ml

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Instructions for use

Intended purpose

For cleansing, moisturing and infection prophylaxis of acute and chronic wounds.

Product features

Wound cleansing and moistening with LAVANOX-Serag® wound spray helps to remove biofilm and fibrinoid debris. This promotes wound healing and reduces the risk of infection. The LAVANOX-Serag® wound spray also supports a low-pain removal of the wound dressing. Wound odours are eliminated by using LAVANOX-Serag® wound spray. The application of LAVANOX-Serag® wound spray is well tolerated and virtually painless. LAVANOX-Serag® wound spray is suitable for wounds on skin and mucous membranes. There are no restrictions on the frequency and duration of use. Residues of LAVANOX-Serag® wound spray may remain in the wound. As a multi-dose container, LAVANOX-Serag® wound spray is preserved with NaOCl/HOCl (<0.08%). This means that the opened single product can be used for three months. A slight chlorine odour is a sign of the activity of the preservative system.

Instructions for use

A wound should first be cleaned of dirt, blood and coatings. For this purpose, the affected areas are wetted with LAVANOX-Serag® wound spray with several sprays from a short distance and, if necessary, cleaned mechanically. If the wound is heavily contaminated, the application can be repeated several times. After the application of LAVANOX-Serag® wound spray, it is not necessary to rinse with another solution (i.e. residues of LAVANOX-Se-

rag® wound spray may remain in the wound). No undesirable side effects are known. Encrusted wound dressings or dressings that are difficult to remove are soaked generously with LAVANOX-Serag® wound spray until they are soaked through. After a short waiting time, these can be carefully removed with minimal pain. After cleaning, a thin layer of wound gel can then be applied if required. The wound gel and the wound spray gel from the LAVANOX-Serag® product line are suitable for this. The frequency of dressing changes depends on the nature of the wound and must be adapted to the respective wound situation. LAVANOX-Serag® wound spray can be used in combination with conventional dressings, plasters and other bandages. The simultaneous use of LAVANOX-Serag® wound spray with pre-soaked wound dressings, antiseptics and other wound care products can weaken their effect. Interactions are possible with chlorhexidine and polyurethane. Dressing materials should therefore be compatible with oxidising agents, especially hypochlorite. However, delayed application is possible. Please also refer to the instructions for use of the respective product.

Patient target group(s)

LAVANOX-Serag® wound spray can be used without restriction with regard to the age and gender of the patient. For use during pregnancy, while breastfeeding and in children under two years of age, consult your doctor or pharmacist beforehand.

Intended users

Health care professionals and laypersons. LAVANOX-Serag® wound spray can be used by medical laypersons if they have been instructed by medical professionals. Children up to the age of 12 should be supervised by an adult.

Medical indications
Acute and chronic wounds.

Limitations of use and contra-indications

LAVANOX-Serag® wound spray should not be used in the case of known allergies to any of the ingredients. If you notice any signs of an allergic reaction to the ingredients, discontinue use of the product. Do not swallow.

Warnings, precautions, measures to be taken


Avoid contact with eyes. In case of eye contact, rinse with plenty of water. Despite the generally very good tolerability, there is a residual risk of a discomforting feeling during application due to the individual wound situation. Allergic reactions cannot be ruled out. Changing the dressing too often (overdosing) can disrupt the wound healing process. The treatment should therefore be monitored by a healthcare professional. If there is redness, swelling or pain in the wound area or if the condition of the wound does not improve within a few days, a doctor should be consulted. Only use the product up to the stated expiry date and within 3 months after opening. Please report any potentially serious incident related to the product immediately to Prontomed (vigilanz@prontomed.de) and the competent authority in your country.


Safe disposal

No special measures are required for safe disposal. Please observe the local disposal guidelines.

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GmbH & Co. KG
Prontomed GmbH