D

# LAVANOX-Serag® Wundsprühgel 75ml

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## Instructions for use

### Intended purpose

For cleansing, long-term moistu-  
rizing and infection prophylaxis  
of acute and chronic wounds.

### Product features

LAVANOX-Serag® wound spray gel  
is a ready-to-use hydrogel for  
painless application to acute and  
chronic wounds. It forms a pro-  
tective film and keeps the wound  
moisturised. It thus forms a  
microenvironment that promo-  
tes wound healing and supports  
the natural healing and cleansing  
process in the wound. The risk  
of infection is reduced. Wound  
odours are eliminated by the  
use of LAVANOX-Serag® wound  
spray gel. Thanks to its proper-  
ties, LAVANOX-Serag® wound  
spray gel can be applied quick-  
ly and easily without the need  
for additional aids and re-  
mains in the wound as a hydro-  
gel. The application of LA-  
VANOX-Serag® wound spray gel  
is well tolerated and virtually  
painless. LAVANOX-Serag® wound  
spray gel is suitable for wounds  
on skin and mucous membranes.  
There are no restrictions on the  
frequency and duration of use.  
Residues of LAVANOX-Serag®  
wound spray gel may remain  
in the wound. As a multi-do-  
se container, LAVANOX-Serag®  
wound spray gel is preserved  
with NaOCl/HOCl (<0.06%). This  
means that the opened single  
product can be used for three  
months. A slight chlorine odour is

of a sign. The activity of the active ingredients should be supervised by an adult. The preservative system.

### Instructions for use

Shake well before use! Before the first application and when changing the dressing, we recommend using a wound spray or wound irrigation solution to remove exudate, fibrinoid deposits and other impurities. Products from LAVANOX-Serag® are suitable for this purpose. Spray LAVANOX-Serag® wound spray gel several times at short intervals until the wound is completely wetted, depending on the condition of the wound. No undesirable side effects are known. The frequency of dressing changes depends on the nature of the wound and must be adapted to the respective wound situation. LAVANOX-Serag® wound spray gel can be used in combination with conventional dressings, plasters and other bandages. The simultaneous use of LAVANOX-Serag® wound spray gel with pre-soaked wound dressings, antiseptics and other wound care products can weaken their effect. Interactions are possible with chlorhexidine and polyurethane. Dressing materials should therefore be compatible with oxidising agents, especially hypochlorite. However, delayed application is possible. Please also refer to the instructions for use of the respective product.

### Patient target group(s)

LAVANOX-Serag® wound spray gel can be used without restriction with regard to the age and gender of the patient. For use during pregnancy, while breastfeeding and in children under two years of age, consult your doctor or pharmacist beforehand.

### Intended users

Health care professionals and laypersons. LAVANOX-Serag® wound spray gel can be used by medical laypersons if they have been instructed by medical professionals. Children up to the age of 12

**Medical indications**  
Acute and chronic wounds.

### Limitations of use and contra-indications

LAVANOX-Serag® wound spray gel should not be used in the case of known allergies to any of the ingredients. If you notice any signs of an allergic reaction to the ingredients, discontinue use of the product. Do not swallow. Do not use gel that has been frozen. Do not use LAVANOX-Serag® wound spray gel on bleeding or heavily oozing wounds.

### Warnings, precautions, measures to be taken


Avoid contact with eyes. In case of eye contact, rinse with plenty of water. Despite the generally very good tolerability, there is a residual risk of a discomforting feeling during application due to the individual wound situation. Allergic reactions cannot be ruled out. Changing the dressing too often (overdosing) can disrupt the wound healing process. The treatment should therefore be monitored by a healthcare professional. If there is redness, swelling or pain in the wound area or if the condition of the wound does not improve within a few days, a doctor should be consulted. Only use the product up to the stated expiry date and within 3 months after opening. Please report any potentially serious incident related to the product immediately to Promotendo (vigilanz@promotendo.de) and the competent authority in your country.


### Safe disposal

No special measures are required for safe disposal. Please observe the local disposal guidelines.

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GmbH & Co. KG
Promotendo GmbH